


## Automated PCR Setup

### Automated Real-Time PCR Set-up on a Tecan Freedom PCR Workstation



#### Introduction

Real-Time chemistry allows for the detection of PCR amplification products during the early phases of the reaction. Measuring the kinetics of the reaction provides a distinct advantage over traditional PCR.

The Tecan Freedom PCR Workstation provides a highly flexible platform for automated set-up of traditional PCR or Real-Time PCR. Using Tecan's disposable tips, precise and accurate set-up even at low volumes can be performed without any risk of cross-contamination.

Moreover, with the integrated robotic manipulator (gripper tool) and appropriate storage modules safe and reliable set-up of multiple batches of 96 samples is possible without any user interaction.


#### PCR Procedure

Multiplex Real-Time PCR can be used for quantification of downy mildew infestation of grape plants. For this purpose, total DNA from plant (*Vitis vinifera*) and parasite (*Plasmopara viticola*) is extracted from infected leaf punches by traditional extraction methods. Multiplex Real-Time

PCR is performed in an ABI PRISM® 7700 instrument (Applied Biosystems) with specific primers for each species, labeled with different fluorescence labels (FAM™ and VIC™, respectively). Typically 36-40 PCR cycles were performed using AmpliTaq® Gold chemistry (Applied Biosystems). To check for potential cross-contamination, negative controls are placed in every second row on the amplification plate.

#### Automation

Automation of the Real-Time PCR assay preparation is realized on a Tecan Freedom PCR Workstation. The instrument is equipped with 8 pipetting channels (disposable tips) and specific carriers for precise positioning of Taqman® plates (Applied Biosystems). Pipetting of mastermix (25 µl) and template DNA (5 µl) is performed with either 200-µl or 50-µl disposable tips (approx. 6 minutes for reaction set-up of 96 samples). Gemini™ 4.0 is used for software control of the process.




**Fig. 1: Assay Set-up**

Full automation of the Real-Time PCR set-up on a Tecan Freedom PCR Workstation using 50  $\mu$ l disposable tips for template pipetting. Multiple plates can be prepared without any user interaction.

To check for potential cross-contamination, a stripe pattern was used for sample (amplification template) and negativ controls (water) distribution into the PCR plate.

## Results

Reliable and precise pipetting is demonstrated by consistent low variation of  $C_t$  values within replicates ( $n=4$ ). All liquid transfers are performed without any cross-contamination as confirmed by non-appearance of amplification products in negative control wells. Reaction set-up of 96 samples and controls is completed within approximately 10 minutes.



**Fig.2: Reliable Real-Time-PCR**

Typical Real-Time amplification profile (Taqman<sup>®</sup>; FAM) of *Vitis vinifera* DNA samples (multiples of three samples (S1, S2, S3) and amplification control with standard template (C) are shown).

$\Delta C_t$  within replicates is consistently  $<1.3$ . No amplification of negative controls was observed.

Analogous results for *Plasmopara viticola* DNA amplification (data not shown).

## Conclusion

The Tecan Freedom PCR Workstation is the ideal instrument for automated Real-Time PCR pipetting:

- high precision pipetting
- no cross-contamination
- deckspace for multiple batches
- flexible access of pipetting tips to various tube and plate types

Complex amplification control patterns are often typical for Real-Time PCR assays. Due to its versatile Liquid Handling Arm (independent z-move-

ment of pipetting tips and flexible tip spacing), the Freedom Sample Preparation Workstation is able to perform any PCR set-up task with high speed and reliability.

Low-volume disposable tips enable precise liquid transfers down to 1  $\mu$ l.

## Acknowledgments

Data kindly provided by Ms. Giorgia Valsesia, Eidgenössische Technische Hochschule, Zurich (Switzerland).

Tecan Asia (Pte) Ltd., 80, Marine Parade, #13-04, Singapore 449269, Singapore, T +65 644 41 886, F +65 644 41 836  
 Tecan Sales Austria GmbH, Untersbergstrasse 1a, A-5082 Grödig / Salzburg, Austria, T +43 62 46 89 33, F +43 62 46 72 770  
 Tecan Sales International GmbH, Untersbergstrasse 1a, A-5082 Grödig / Salzburg, Austria, T +43 62 46 89 33, F +43 62 46 72 770  
 Tecan Benelux B.V.B.A., Vaartdijk 55, B-2800 Mechelen, Belgium, T +32 15 42 13 19, F +32 15 42 16 12  
 Tecan Benelux B.V.B.A., Industrieweg 30, NL-4283 GZ Giessen, Netherlands, T +31 18 34 48 17 4, F +31 18 34 48 06 7  
 Tecan Deutschland GmbH, Theodor-Storm-Straße 17, D-74564 Crailsheim, Germany, T +49 79 51 94 170, F +49 79 51 50 38  
 Tecan France S.A., Parc d'Activités de Pissaloup, Bâtiment Hermes II, Rue Edouard Branly, F-78190 Trappes, France, T +33 1 30 68 81 50, F +33 1 30 68 98 13  
 Tecan Italia S.r.l., Via F.lli Cervi, Palazzo Bernini, Centro Direzionale Milano 2, I-20090 Segrate (Mi), Italy, T +39 02 215 21 28, F +39 02 215 97 441  
 Tecan Japan Co. Ltd, Meiji Seimei Fuchu Building 10F, 1-40 Miyamachi, Fuchu City, Tokyo, Japan, T +81 42 334 88 55, F +81 42 334 04 01  
 Tecan Nordic AB, Box 208, SE-431 23 Mölndal, Sweden, T +46 31 75 44 000, F +46 31 75 44 010  
 Tecan Portugal, Quinta da Fonte - Edifício Pedro I, 2780-730 Paço d'Arcos, Portugal, T +351 21 000 82 16, F +351 21 000 16 75  
 Tecan Sales Switzerland AG, Seestrasse 103, CH-8708 Männedorf, Switzerland, T +41 1 922 89 22, F +41 1 922 89 23  
 Tecan Spain, Sabino de Arana, 32, E-08028 Barcelona, Spain, T +34 93 490 01 74, F +34 93 411 24 07  
 Tecan UK, Theale Court, 11-13 High Street, Theale, UK-Reading RG7 5AH, United Kingdom, T +44 11 89 300 300, F +44 11 89 305 671  
 Tecan US, P.O. Box 13953, Research Triangle Park, NC 27709, USA, T +1 919 361 5200, F +1 919 361 5201